


BASIC INFORMATION

Straumann® RidgeFit Implants System



ABOUT THIS GUIDE

This surgical and prosthetic procedure describes the steps required for implantation and restoration of the Straumann® RidgeFit Implants System. The Straumann® RidgeFit Implants System is recommended for use only by clinicians with advanced surgical skills. It is assumed that the user is familiar with placing dental implants. Not all detailed information will be found in this guide. Reference to existing Straumann procedure manuals will be made throughout this document.

Not all products shown are available in all markets.

CONTENTS

1. THE STRAUMANN® RIDGEFIT IMPLANTS SYSTEM	3
1.1 Portfolio overview	4
1.2 The Straumann® RidgeFit Implants at a glance	5
1.3 Straumann® Optiloc® Retentive System	6
2. SURGICAL PROCEDURE	7
2.1 Preoperative planning	7
2.2 Instruments	9
2.3 Implant bed preparation	13
2.4 Implant insertion	14
3. PROSTHETIC PROCEDURE	18
3.1 Chairside modification of an existing well-fitting and well-functioning lower denture into an overdenture supported by Optiloc® Retentive System/Straumann® RidgeFit Implants	18
3.2 Creating a new overdenture with the Optiloc® Retentive System	21
4. USING THE OPTILOC® TOOLS	24
4.1 Optiloc® Matrix Housing Extractor	24
4.2 Optiloc® Laboratory Instrument (blue)	24
4.3 Optiloc® Retention Insert Instrument	25
5. SPECIALLY FEATURED OPTILOC® COMPONENTS	26
6. PRODUCT REFERENCE LIST	27
6.1 Straumann® RidgeFit Implants Roxolid® SLA®	27
6.2 Optiloc® Processing Package, Retention Inserts and Matrix Housings	27
6.3 Optiloc® tools and auxiliary parts	28
6.4 Straumann® Modular Cassette	28
7. FURTHER INFORMATION	29

1. THE STRAUMANN® RIDGEFIT IMPLANTS SYSTEM

The Straumann® RidgeFit Implants System offers one-piece Tissue Level implants with an Optiloc® prosthetic connection. These are designed for narrow edentulous ridges and immediate treatment procedures (if at least 35 Ncm insertion torque is achieved in all implants) to stabilize full removable overdentures.

The Straumann® RidgeFit Implants are made from the material Roxolid® with the SLA® surface and are available in the endosteal diameters Ø 2.4 mm, with length options of 10 mm, 12 mm and 14 mm.






To obtain more information about indications and contraindications related to the implant, please refer to the corresponding instructions for use. Instructions for use can be found on www.ifu.straumann.com

For further information on the Optiloc® Retentive System please refer to www.ifu.valoc.ch

1.1 PORTFOLIO OVERVIEW

Surgical components

					
RidgeFit Implants GH 2.8 mm: 042.947S, 10 mm 042.948S, 12 mm 042.949S, 14 mm GH 3.8 mm: 042.954S, 10 mm 042.955S, 12 mm GH 4.8 mm: 042.956S, 10 mm 042.957S, 12 mm	Needle Drill, long 027.0007S	Pilot Drill 027.0011S	Adapter Optiloc® for Ratchet 170.2	Adapter Optiloc® for Handpiece 170.1	Paralleling Posts 046.796

Prosthetic components

				
Optiloc® Model Analogue 2102.0024-OPT	Matrix Housings 2102.0001-OPT 2102.0009-OPT 2102.0010-OPT	Retention Insert 2102.0003-OPT 2102.0004-OPT 2102.0005-OPT 2102.0006-OPT 2102.0007-OPT 2102.0008-OPT	Optiloc® Block Out Spacer 2102.0023-OPT	Optiloc® Processing Collar 2102.0011-OPT
				
Optiloc® Impression Coping 2102.0012-OPT	Optiloc® Matrix Housing Extraction Instrument 3202.0003-OPT	Optiloc® Laboratory Instrument 3202.0002-OPT	Optiloc® Retention Insert Instrument 3202.0001-OPT	Equipment Box with 3 Instruments 5102.0000-OPT

1.2 THE STRAUMANN® RIDGEFIT IMPLANTS AT A GLANCE

Optiloc®

- Minimized maintenance, narrow diameter
- TiN coating for good wear resistance

Apically tapered implant body design allows underpreparation and supports a high primary stability



Roxolid®:

- High material strength and biocompatibility
- Peace of mind with RidgeFit Implants


SLA®:

- Predictability in osseointegration
- Scientific evidence
- Low prevalence of peri-implantitis
- Bone preservation

1.2.1 Material

Roxolid® is a groundbreaking material specifically designed for the use in dental implantology. The titanium-zirconium alloy is stronger than pure titanium^{1,2} and has excellent osseointegration properties³⁻⁵. This combination of properties is unique in the market, there is no other metallic alloy which unifies high mechanical strength and osteoconductivity.

Thanks to their outstanding biological and mechanical properties, Roxolid® Implants offer more treatment options than conventional titanium implants.



Roxolid® shows a 20% higher tensile strength than Straumann cold worked titanium and a 80% higher strength than standard titanium Grade 4.

1.3 STRAUMANN® OPTILOC® RETENTIVE SYSTEM

The Straumann® Optiloc® Retentive System for full removable overdentures offers an innovative connection coating (TiN) with an excellent wear resistance, overcoming up to 40° implant convergence or divergence. Together with its durable PEEK¹ matrices the Optiloc® Retentive System provides a unique and long-lasting attachment performance.

1.3.1 Straumann® Optiloc® Retentive System at a glance

1


- PEEK¹ Matrix Housings offering excellent chemical and physical properties
- Matrix accommodates up to 40° prosthetic divergence between two abutments
- 6 retention strengths offer optimal adjustment of the denture retention
- Matrix Housing available in titanium, or color-neutral PEEK¹ for a more aesthetic outcome

2

- Offering a smooth surface and ultimate hardness for excellent wear resistance



The Optiloc® Matrix System allows a convergence, or divergence, of up to 20 degrees of each implant in relation to the denture's path of insertion. This means that divergences between two implants of a maximum of 40 degrees can be corrected.



¹ Polyether ether ketone

2. SURGICAL PROCEDURE

The workflow for the surgical procedure for the Straumann® RidgeFit Implants System involves 3 steps:


- Preoperative planning
- Implant bed preparation and
- Implant insertion.

2.1 PREOPERATIVE PLANNING

After patient selection and evaluation protocols have been completed, the number of Straumann® RidgeFit Implants that should be placed (minimum of four in the mandible, minimum of six in the maxilla) are determined and thoroughly discussed with the patient. Information on bone availability for the implant bed of the patient and information of tissue depth mucosa thickness in the region of the prospected implant site by measuring with a perio probe should be available. After site selection, Straumann® RidgeFit Implants should be placed at least 5 mm apart.

When anatomic conditions allow, distribute the number of implants along the arch in-order to minimize cantilever and to provide optimal load distribution/better load conditions. When anatomic situation is not optimal for the mandible the implants should be placed beginning at least 5 mm anterior to the mental foramen. The remaining anterior space should be distributed equally between implants and respecting the minimum distance between implants (5mm).


For maxillary placement, careful implant length selection must be followed to avoid anatomical structures such as nasal cavity and maxillary sinus.



2.1.1 X-ray reference foil

The vertical bone availability determines the maximum allowable length of the implant that can be placed. For easier determination of the vertical bone availability, we recommend the use of an x-ray reference foil with X-ray Reference Sphere (Art. No. 049.076V4).

Similar to the distortions that occur in X-rays, the implant dimensions are shown on the individual templates with the corresponding distortion factors (1:1 to 1.7:1). Determining each magnification factor or scale is facilitated by showing the X-ray reference sphere on the template (next to the scale reference).



Note: Use only the x-ray template specific to the implant type.

To calculate the effective bone availability, use the following formula:

$$\frac{\text{X-ray Reference sphere 5 mm} \times \text{bone availability (X-ray*)}}{\text{Reference sphere diameter on the X-ray}} = \text{effective bone availability}$$

2.1.2 Surgical preparation

	Bone Type				Soft tissue depth		Buccolingual width		
	Type I	Type II	Type III	Type IV	< 2 mm	≥ 2 mm	< 4.4 mm	≥ 4.4 mm with flap	≥ 5.4 mm flapless**
2.4 mm Straumann® RidgeFit Implants	✓	✓	✓	✗	✗	✓	✗	✓	✓

✗ Not recommended ✓ Recommended

* Taking into consideration all implant-related anatomic structures (e.g. mandibular canal, sinus maxillaris, etc.)

** Flapless procedures have a higher planning inaccuracy. We recommend at least a ridge width of 5.4 mm for such interventions

This implant is contraindicated for bone class 4. Additionally it is recommend that the treatment should only be conducted with patients that have more then 2mm soft tissue depth and a buccolingual width of more than 4.4 mm. The Procedure can be done flap-less only if there is more then 5.4mm buccolingual bone width.

2.2 INSTRUMENTS

2.2.1 Drills

The Straumann® instruments have depth marks at 2 mm intervals that correspond to the available implant lengths. The first bold mark on the drills represents 10 mm and 12 mm, where the lower edge of the mark corresponds to 10 mm and the upper edge to 12 mm. The drills are delivered sterile.



1. Drill: 027.00075
2. Drill: 027.00115
3. Implant: 042.9475


2.2.2 Paralleling Post

The Paralleling Post is an instrument used to ensure the correct, parallel positioning of the implant during implant bed preparation and to align with other implants.

As a secondary feature, the mid portion of the Paralleling Post represents the gingiva height/machined part of the implant.


Characteristics

- Ø 1.6 mm
- Ø 2.2 mm
- Height 2.8 mm (gingiva height implant)
- Material: TAV
- Delivered sterile






2.2.3 Adapter

Specific adapter to use for insertion of the Straumann® RidgeFit Implants.



2.2.4 Ratchet and Torque Control Device

The Ratchet is a two-part lever arm instrument with a rotary knob for changing the direction of force. It is supplied with a service instrument, which is used to tighten and loosen the head screw. The Holding Key (046.064) can be used to stabilize the Ratchet.


Ratchet and Torque Control Devices			
	Holding Key	Ratchet	BLX Torque Control Device for Ratchet, Surgical
			
Intended use	Auxiliary	Torque transmission	Surgical
Torque markings	NA	NA	0 / 35 / 50 / 80 Ncm
Article Number	046.064	046.119	066.1100
Material	Stainless steel	Stainless steel	Stainless steel, DLC coated

Note: To ensure prolonged perfect function, the Ratchet must always be taken apart and the individual parts disinfected, cleaned and sterilized after use. Its function must be checked in good time before each use.

2.2.5 Setup for Straumann® Mini Implant freehand surgery

For more information refer to *Straumann® Modular Cassette Selection Guide* (702824/en).


A Module
041.761



Caution: these instruments are for single use only and can not be sterilize with the Modular Cassette and other instruments.

2.3 IMPLANT BED PREPARATION

2.3.1 Drilling protocol for Straumann® RidgeFit Implants



Recommended speed: rpm max 800

2.4 IMPLANT INSERTION



Stabilization of a mandibular denture

A minimum of four Straumann® RidgeFit Implants should be placed in the mandible.

Caution: Pay attention to the Inferior Alveolar Nerve, and the sub-lingual artery.

Note: Always start with the most distal implant at least 7 mm anterior to the mental foramen.



Step 1 – Site preparation (flapless)

Entry point is marked on the patient's tissue (tissue punch optional whenever there is sufficient attached mucosa available). No incision is necessary with this procedure. This punching procedure is only recommended if sufficient attached mucosa remains around the implant to ensure peri-implant health long term.

Mark the implantation site determined during the implant position planning with the Ø1.6 mm Needle Drill.


In presence of a thin ridge, the use of a round bur may be necessary and shall be used in order to mark the bone before using the 1.6 mm Drill.



Step 2 – Implant axis

Mark the implant axis with the Needle Drill to a depth of 6 mm.


For Ø 2.4 mm RidgeFit Implants in **soft bone (type III)**, the implant bed preparation ends here.



Drill the implant bed to the final depth with the 1.6 mm Needle Drill, while correcting unsatisfactory implant axis orientation if necessary. Use the 1.6 mm side of the Paralleling Post to check the implant axis.

For \varnothing 2.4 mm RidgeFit Implants in **hard bone (type II)**, the implant bed preparation ends here. (Optionally can be continued with the 2.2mm BLT Pilot Drill).

Note: In case of vertically reduced bone availability, an x-ray should be taken at this step in order to be sure that drill did not pass through the mandibular basal bone.



Step 3 – Optional: Widen implant bed to \varnothing 2.2mm in type II bone
With the \varnothing 2.2mm BLT pilot Drill, drill to a depth of about 6mm.

In very hard bone (type I): Drill the implant bed to the final depth with the 2.2mm BLT pilot drill.




Insert the \varnothing 2.2mm Paralleling post to check for correct implant axis orientation. Use the \varnothing 2.2mm BLT Pilot Drill to prepare the implant bed to final preparation depth.



Further implants, implant alignment:

Leave in the paralleling post and proceed with the next implant bed preparation and continue until all implants are placed. Always orienting the drill to the most proximal implant bed with the paralleling post. Distribute the implants equally, respecting the minimum distance between implants (5 mm).

Note: Parallelism of the implants is essential.



Step 4 – Implant placement

Straumann® RidgeFit Implants are delivered in a sterile vial and mounted on the vial cap which serves as initial insertion tool.

Note: When opening the sterile blister, ensure that the sterile vial is securely seated on a sterile surgical tray.



Use of the vial cap as finger driver

Insertion of implant begins with the vial cap until more torque is necessary.

Caution: Please do not use implant if it is detached from the vial cap after opening of the blister.



Step 5 – Final implant positioning

Place implant

Straumann® RidgeFit Implants can be placed with the Handpiece or manually with the Ratchet. A maximum speed of 15 rpm is recommended.

Use the Ratchet and/or Handpiece to move the implant into its final position turning it clockwise.



Final placement is achieved once the entire conditioned SLA® surface is engaged into the bone

Note: For immediate loading a minimum insertion torque of 35Ncm is necessary for all implants. If 35Ncm is not achieved on all implants conventional loading is recommended. Early loading is contraindicated in all cases. For further information please go to prosthetic workflow.

Do not exceed 80 Ncm insertion torque during implant placement as this may lead to implant damage.

Note: A minimum of 4 Straumann® RidgeFit Implants should be placed to stabilize a full mandibular denture and a minimum of 6 Straumann® RidgeFit Implants should be placed to stabilize a full maxillary denture.



Maxillary denture stabilization

Proceed as above but pay special attention to:

A minimum of six Straumann® RidgeFit Implants should be placed in the maxilla

3. PROSTHETIC PROCEDURE

3.1 CHAIRSIDE MODIFICATION OF AN EXISTING WELL-FITTING AND WELL-FUNCTIONING LOWER DENTURE INTO AN OVERDENTURE SUPPORTED BY OPTILOC® RETENTIVE SYSTEM/STRAUMANN® RIDGEFIT IMPLANTS

Caution: It is a prerequisite, however, that the lower complete denture does not need to be relined by a dental technician.

This work-flow illustration refers to a lower denture procedure. The upper denture follows the same steps.



Step 1 – Place white processing collar on each Optiloc®

The Processing Collar is used to block out the area surrounding the Optiloc®.

Then place a Matrix Housing with a Retention Insert (recommendation white, medium) onto each Optiloc® abutment, leaving the white processing collar beneath it.



Step 2 – Prepare the complete denture

Hollow out the existing denture base in the areas of the Optiloc® Matrix Housings with Handpiece and resin bur. There should be a minimum space of 1 mm around the housings to allow for sufficient thickness of the self-polymerizing resin.



Step 3 – Seat denture

Use wash impression silicone to confirm adequate clearance between the Matrix Housings and the denture base.

Insert the complete denture into the patient's mouth and check the clearance. The Matrix Housings fixed on the abutments should not touch the denture base. Reconfirm adequate space using wash impression silicone. Adjust the denture base until seated passively in occlusion without touching the Matrix Housings.



Step 4 – Prepare denture

Apply proper adhesive/primer material in the denture prior to receiving the self-curing PMMA resin.



Step 5 – Polymerize the Matrix Housings

Fill the hollowed area with self-curing PMMA resin to polymerize the Matrix Housings in the denture.

Apply a small amount of acrylic resin to the tissue-contact surface of the denture and around the Matrix Housings. Insert the complete denture into the oral cavity.



Step 6 – Seat denture in occlusion

Once the complete denture is properly seated, maintain the patient in centric occlusion while the acrylic sets.



Step 7 – Discard Optiloc® Processing Collar

Once the resin has cured, remove the complete denture from the mouth and discard the white Optiloc® Processing Collar.

Put the complete denture in hot, but not boiling, water. Place it in a pressure pot when available.



Step 8 – Finish denture

After final curing, remove any excess acrylic and finish the denture base.

If needed, exchange the white, medium Optiloc® Retention Insert with other Optiloc® Retention Inserts and insert the final overdenture into the patient's mouth.

3.2 CREATING A NEW OVERDENTURE WITH THE OPTILOC® RETENTIVE SYSTEM

Procedure in the dental office – Impression taking on abutment level

This work-flow illustration refers to a lower denture procedure. The upper denture follows the same steps.




Step 1 – Placing the Optiloc® Impression Coping
Place the Impression Coping on the Optiloc®



Step 2 – Impression taking
Use the mucodynamic technique for impression taking (vinyl polysiloxane or polyether rubber).

Send the impression to the dental lab.


Procedure in the dental lab



Step 1 – Inserting the Optiloc® Model Analogue
Insert the Optiloc® Model Analogue into the Optiloc® Impression Coping (see chapter 4 using the Optiloc® tools)



Step 2 – Fabricating the master cast
Pour a master model using standard methods and type-4-dental stone (DIN 6873).



Step 3 – Placing the Optiloc® Processing Collar and Matrix Housing
Place the Matrix Housing incl. a Retention Insert (e.g. 2102.0005-STM, white, medium) onto the Optiloc®.

For a chairside polymerization of the Optiloc® Matrix Housing use the Optiloc® Block Out Spacer to create the space needed.



Step 3.1 – Finalizing the new Optiloc overdenture

Place white Processing Collar on all Optiloc® Model Analogues.



Step 3.2 – Processing the overdenture

Process the overdenture according to the standard procedures.



The dental lab will return the finalized Optiloc® overdenture to the dental office.

Procedure in the dental office



Step 4 – Seating the new Optiloc® overdenture

Select the appropriate Optiloc® Retention Insert (see chapter 5 Special featured Optiloc® components).



Step 4.1 – Selecting and inserting the Optiloc® Retention Inserts

Exchange the Optiloc® Retention Inserts to the Matrix Housing using the Retention Insert Instrument (brown) (see chapter 4 Using the Optiloc® Tools).



Step 4.2 – Seating the finished overdenture

Seat the finished overdenture.

Soft Reline Protocol when primary stability is not archived on all implants.

1. Grind down denture base from the existing denture. At least 1mm and relieve denture to accommodate the prosthetic heads of each implant. Important: At this step the implants heads must be absent of any contact.
2. Roughen and degrease the tissue-contact surface appropriately.
3. Apply soft relining material onto the tissue-contact surface of the denture.
4. Place the denture on the patient's mouth and ask patient to apply normal bite pressure in centric occlusion.
5. Allow proper setting time according to the relining material brand of choice.
6. Remove denture and trim excess material with fine scissors or a surgical blade. When available apply glazing material.
7. Do not remove the palate of a maxillary denture during this stage.
8. Ask Patient to keep the denture in place for the first 48 hours after placement to prevent tissue overgrowth.
9. With a healing time minimum of two months the soft reliner material is replaced with the final prosthesis. Please refer to Prosthetic workflow 3.1 and 3.2
8. After osseointegration the palatal plate in maxillary denture can be progressively removed, if desired

4. USING THE OPTILOC® TOOLS

4.1 OPTILOC® MATRIX HOUSING EXTRACTOR (FIG. 1)


Removing the Optiloc® Matrix Housing from an overdenture

1. Heat the Optiloc® Matrix Housing Extraction Instrument head (Fig. 2).
2. Apply the hot Optiloc® Matrix Housing Extraction Instrument to the Matrix Housing and let the heat transfer for 2–3 seconds melting the resin around the Matrix Housing.
3. Tilt the Optiloc® Matrix Housing Extraction Instrument to the opposite side of the beak-shape end to remove the Optiloc® Matrix Housing (Fig. 3).

4.2 OPTILOC® LABORATORY INSTRUMENT (BLUE) (FIG. 4)

Placing the Optiloc® Model Analogue

1. Pick up the Optiloc® Model Analogue with the opposite side of the Optiloc® Laboratory Instrument (Fig. 5/6).
2. Position the Optiloc® Model Analogue in the impression (Fig. 7).





4.3 OPTILOC® RETENTION INSERT INSTRUMENT (FIG. 8)

Mounting the Optiloc® Retention Insert

1. Pick up the Optiloc® Retention Insert with the gripper end of the Optiloc® Retention Insert Instrument. The Optiloc® Retention Insert will lock on to the tool (Fig. 9).
2. Place the Optiloc® Retention Insert into the Optiloc® Matrix Housing (Fig. 10). The Optiloc® Retention Insert “clicks” into position (Fig. 11).

Demounting the Optiloc® Retention Insert

1. Apply the plunger end of the Optiloc® Retention Insert Instrument to the Optiloc® Retention Insert and engage with light pressure (Fig. 12/13).
2. Remove the Optiloc® Retention Insert from the Optiloc® Matrix Housing using a slight rotational movement (Fig. 14).
3. Use the special indentation in the handle of the Optiloc® Matrix Housing Extraction Instrument (Fig. 1) to remove the Optiloc® Retention Insert from the Optiloc® Retention Insert Instrument with a tilting movement (Fig. 15/16).



5. SPECIALLY FEATURED OPTILOC® COMPONENTS



Optiloc® Retention Inserts

The Optiloc® Matrix System allows a convergence, or divergence, of up to 20 degrees of each implant in relation to the denture's path of insertion.


Note:

It is recommended to use the light retention force first (white). In case it feels too loose for the patient, exchange with inserts with a higher retention force.



Optiloc® Processing Collar

The Processing Collar blocks out the area surrounding the abutment, preventing resin or a bonding agent from flowing into the Matrix Housing and embedding the abutment.



Optiloc® Matrix Housing - Extended

This Matrix Housing - Extended offers an extended attachment option. It is used for low-lying abutment heights or in situations requiring more retention. The attachment may be shortened according to the required height.

















Optiloc® Block Out Spacer











The Optiloc® Block Out Spacer is a placeholder for the Optiloc® Matrix Housing. It is used for the model cast, cast metal-reinforced denture or if the Optiloc® Matrix Housing is to be polymerized into the overdenture chairside.

6. PRODUCT REFERENCE LIST













6.1 STRAUMANN® RIDGEFIT IMPLANTS ROXOLID® SLA®

Art. No.		Article
Straumann® RidgeFit Implants		
042.947S		RidgeFit Implants Ø 2.4 mm, GH 2.8 mm, SLA®, TiN, 10 mm, Roxolid®
042.954S		RidgeFit Implants Ø 2.4 mm, GH 3.8 mm, SLA®, 10 mm, TiN, Roxolid®
042.956S		RidgeFit Implants Ø 2.4 mm, GH 4.8 mm, SLA®, 10 mm, TiN, Roxolid®
042.948S		RidgeFit Implants Ø 2.4 mm, GH 2.8 mm, SLA®, TiN, 12 mm, Roxolid®
042.955S		RidgeFit Implants Ø 2.4 mm, GH 3.8 mm, SLA®, 12 mm, TiN, Roxolid®
042.957S		RidgeFit Implants Ø 2.4 mm, GH 4.8 mm, SLA®, 12 mm, TiN, Roxolid®
042.949S		RidgeFit Implants Ø 2.4 mm, GH 2.8 mm, SLA®, TiN, 14 mm, Roxolid®
Auxiliary Parts		
046.796		Paralleling Post for Mini Implants, sterile
170.1/170.1-OPT		Adapter Optiloc® for Handpiece, length 26 mm
170.2/170.2-OPT		Adapter Optiloc® for Ratchet, length 17 mm
027.0007S		Needle Drill, long, single use
027.0011S		BLT Pilot Drill, long, Ø 2.2 mm, single-use, TAN
2102.0024-STM/ 2102.0024-OPT		Optiloc® Model Analogue, blue, 4 pcs.
2102.0012-STM/ 2102.0012-OPT		Optiloc® Impression Coping, white, 4 pcs.





6.2 OPTILOC® PROCESSING PACKAGE, RETENTION INSERTS AND MATRIX HOUSINGS

Art. No.		Article
Processing Package		
5202.0001-STM/ 5202.0001-OPT		Optiloc® Processing Package Optiloc® Matrix Housing, titanium, 2 pcs. Optiloc® Retention Insert, white, light, 2 pcs. Optiloc® Retention Insert, yellow, medium, 2 pcs. Optiloc® Retention Insert, green, strong, 2 pcs. Optiloc® Processing Collar, silicone, 2 pcs.
Retention Inserts		
2102.0003-STM/ 2102.0003-OPT		Optiloc® Retention Insert, red, extra-light, 4 pcs.
2102.0004-STM/ 2102.0004-OPT		Optiloc® Retention Insert, white, light, 4 pcs.
2102.0005-STM/ 2102.0005-OPT		Optiloc® Retention Insert, yellow, medium, 4 pcs.
2102.0006-STM/ 2102.0006-OPT		Optiloc® Retention Insert, green, strong, 4 pcs.
2102.0007-STM/ 2102.0007-OPT		Optiloc® Retention Insert, blue, extra-strong, 4 pcs.
2102.0008-STM/ 2102.0008-OPT		Optiloc® Retention Insert, black, ultra-strong, 4 pcs.
Matrix Housings		
2102.0001-STM/ 2102.0001-OPT		Optiloc® Matrix Housing, 4 pcs.
2102.0009-STM/ 2102.0009-OPT		Optiloc® Matrix Housing - Elliptic, 4 pcs.
2102.0010-STM/ 2102.0010-OPT		Optiloc® Matrix Housing - Extended, 4 pcs.

6.3 OPTILOC® TOOLS AND AUXILIARY PARTS

Art. No.		Article
5102.0000-STM/ 5102.0000-OPT		Optiloc® Equipment Box, with 3 Instruments Optiloc® Laboratory Instrument (blue) Optiloc® Retention Insert Instrument (brown) Optiloc® Matrix Housing Extractor (gray)
2102.0023-STM/ 2102.0023-OPT		Optiloc® Block Out Spacer, white, 4 pcs.
2102.0011-STM/ 2102.0011-OPT		Optiloc® Processing Collar, silicone, 10 pcs.
3202.0001-STM/ 3202.0001-OPT		Optiloc® Retention Insert Instrument
3202.0002-STM/ 3202.0002-OPT		Optiloc® Laboratory Instrument
3202.0003-STM/ 3202.0003-OPT		Optiloc® Matrix Housing Extraction Instrument
046.795		X-ray Reference Foil for RidgeFit Implants
049.076V4		X-ray Reference Spheres, Ø 5 mm, stainless steel
046.119		Ratchet, including service instrument, length 84 mm, stainless steel
066.1100		Torque Control Device for Ratchet – surgical, stainless steel
046.064		Holding Key, length 85 mm, stainless steel
045.111V4		Cleaning Brush for Ratchet, length 100 mm, Ø 4.5 mm, stainless steel/nylon

6.4 STRAUMANN® MODULAR CASSETTE

Art. No.		Article
041.761		Straumann® Modular Cassette, A Module
041.766		A Module, Ratchet Tray
041.764		A Module, Grommet Tray, 3 small + 3 large
041.762		A Module, Grommet Tray, 6 small

7. FURTHER INFORMATION

For more detailed information on the instructions for use, please consult the following documents:

- *Instructions for Use: Straumann® RidgeFit Implants Ø 2.4 mm Roxolid® SLA® TiN* (704053) <http://ifu.straumann.com>
- *Optiloc® Instructions for Use* <http://ifu.valoc.ch/>
- *Straumann® Surgical and Prosthetic Instruments, Care and Maintenance* (702000/en)
- *Instructions for Use: Straumann® Novaloc® and Optiloc® Martrix System* (704486) <http://ifu.straumann.com>
- *Straumann® Modular Cassette, Basic Information* (702527/en)
- *Instructions for Use: Straumann® Non-sterile Surgical Instruments and Auxiliaries* (701124) <http://ifu.straumann.com>
- *Instructions for Use: Straumann® Drills and Instruments* (701124) <http://ifu.straumann.com>
- *Instructions for Use: Straumann® Prosthetic Planning and Placement Tools* (702879) <http://ifu.straumann.com>
- *Instructions for Use: Straumann® Impression Components* (703287) <http://ifu.straumann.com>

¹ Norm ASTM F67 (states min. tensile strength of annealed titanium). ² Data on file for Straumann cold-worked titanium and Roxolid® Implants, MAT 13336, 20131009. ³ Gottlow J et al. : Evaluation of a new titanium-zirconium dental implant: a biomechanical and histological comparative study in the mini pig. *Journal of Clinical Implant Dentistry and Related Research* 2012; 14: 538-545 ⁴ Wen B et al. : The osseointegration behavior of titanium-zirconium implants in ovariectomized rabbits. *Clin Oral Implants Res.* 2013 Feb 21. ⁵ Barter S et al. : A pilot study to evaluate the success and survival rate of titanium-zirconium implants in partially edentulous patients: results after 24 months of follow-up. *Clin Oral Implants Res.* 2012 Jul;23(7):873-81

International Headquarters

Institut Straumann AG

Peter Merian-Weg 12

CH-4002 Basel, Switzerland

Phone +41 (0)61 965 11 11

Fax +41 (0)61 965 11 01

www.straumann.com

© Institut Straumann AG, 2023. All rights reserved.

Straumann® and/or other trademarks and logos from Straumann® mentioned herein are the trademarks or registered trademarks of Straumann Holding AG and/or its affiliates.

