ES

FILLAH

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FILLAH

CROSSLINKED HYALURONIC ACID
INJECTABLE HYDROGEL

DESCRIPTION
FILLAH is a biocompatible hydrogel of crosslinked sodium hyaluronate obtained by biofermentation and highly purified, dissolved in a physiological and isotonic buffer. The product is an injectable hydrogel, sterile and apyrogenic, with a low level of endotoxins.

COMPOSITION
FILLAH is formulated with highly purified sodium hyaluronate, pharmacopeia quality salts and water for injection (WFI):

MODELS	FILLAH SOFT 1.5%	FILLAH MEDIUM1.75%	FILLAH INTENSE 2%
Composition	(mg/mL)	(mg/mL)	(mg/mL)
Sodium hyaluronate cross-linked	15	17.5	20
Disodium phosphate dodecahydrate	0.6	0.6	0.6
Sodium dihydrogen phosphate dihydrate	0.05	0.05	0.05
Sodium chloride	8	8	8
Water For Injection (WFI)	q.s.	q.s.	q.s.

PACKAGING
FILLAH is supplied in a graduated, pre-filled and disposable Type I glass sterile syringe with a Luer-lock system, containing 1 mL of solution. The product includes as an accessory 2 needles of 27G and 13mm. The manufacturer recommends the administration of the product with the supplied needle or one with identical characteristics.

INTENDED USE
FILLAH is intended as a dermatological treatment for aesthetic purposes. Administered by intradermal injection, the product is indicated to restore and correct facial volume and face contouring and fill in wrinkles and folds.

PROPERTIES
FILLAH restores the loss of volume and face contouring. FILLAH fills in the folds and wrinkles that appear with age. FILLAH redefines the contour and brings fullness to the lips.

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DOSAGE AND ADMINISTRATION
FILLAH must be administered by intradermal injection by qualified medical practitioner following all the rules regarding the aseptic procedure and the injection technique. FILLAH should be injected with a suitable sterile needle (27G and 13mm needles are recommended). It is recommended to administer applying the lowest necessary pressure and massage the treated area after the injection in order to guarantee the uniform distribution of the product. The amount of product to be injected as well as the selection of the most suitable model, will be decided by the specialist based on the area to be treated and the desired level of correction. Any unused portion of the syringe should be discarded.

- DIRECTIONS FOR USE**
- Prior to the injection, the area to be treated should be disinfected thoroughly and bring the product to room temperature.
 - Peel off the sterile barrier and remove the syringe from inside.
 - Remove the tip cap of the syringe without touching the syringe tip in order to maintain it sterile.
 - Attach a suitable sterile needle (27G recommended) to the tip of the syringe. Verify that the needle is properly locked to the Luer-lock adapter of the syringe before using it.
 - Remove air from the syringe, if necessary.
 - Inject the product by intradermal injection into the area to be treated.
 - After finishing the treatment session, discard the syringe and the remaining product.

CONTRAINDICATIONS
Do not administer FILLAH to patients with known hypersensitivity (allergy) to hyaluronic acid and / or other components of the product.

Do not administer FILLAH to patients with acute or chronic dermatological disorders or an inflammatory and / or infectious type skin problem.

Do not administer FILLAH to patients with coagulation defects and / or who are treated with anticoagulant medication.

Do not administer FILLAH in immediate association with laser or pulsed light treatments.

FILLAH should not be administered to pregnant or breastfeeding women.

FILLAH should not be administered to children.

WARNINGS AND PRECAUTIONS
FILLAH should not be injected into blood vessels or eyelids; intravascular injection may produce embolization, occlusion of the vessels, ischemia, blindness or infarction. FILLAH should be administered exclusively by intradermal injection following a strict aseptic procedure.

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CE 0318 

FILLAH should be administered with sterile Luer-lock needles (as supplied with the product). It is recommended that the patient do not wear makeup for 12 hours after the injection and avoid prolonged exposure to the sun, to ultraviolet rays, at temperatures below 0°C, as well as the sauna or Turkish bath during two weeks following the injection. The product should not be used if the package is damaged or broken. The product should not be used after the expiration date. The product is for single use. No reuse; risk of infection. Do not re-sterilize, re-sterilization can affect the physicochemical properties of the product and compromise its effectiveness and safety. The syringe and the used needle should be discarded in a container provided for this purpose. Keep out of children's reach.

INTERACTIONS
Do not use with other products containing quaternary ammonium salts as benzalkonium chloride or with any surgical material treated with this type of substances because sodium hyaluronate can precipitate in their presence. There is no known interaction with other local anesthetics.

ADVERSE EVENTS
The most common side effects are: bruising, nodules or lumps at the injection site, swelling, pain, oedema and redness of the injected area. These effects usually disappear by themselves a few days after the injection. To a lesser extent, cases of immediate or delayed hypersensitivity have been described after injections with hyaluronic acid.

STORAGE
Protect from light and avoid freezing. Store between 5°C and 25°C.

SELL BY DATE
Do not use after the date marked on the package.

MANUFACTURED BY
 **i+Med S. Coop.**
Parque Tecnológico de Álava
Calle de Albert Einstein, 15. Pabellón 15
01510 Vitoria-Gasteiz, (Álava - Spain)
www.imasmed.com

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Segregation and waste management according to current legislation

CE 0318 

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**LABORATORIOS NORMON, S.A.**